

FITTING GUIDE

- STEP 1:** Update spectacle refraction and Add power
- STEP 2:** Select **toric lens** power (adjusted for vertex distance if necessary)
- STEP 3:** Select **Add power** following guidance below:

ADD SELECTION:

SPECTACLE ADD	BOTH EYES
+0.75D to +1.50D	Low Add
+1.75D to +2.50D	High Add

See package insert for more information.

EVALUATE THE LENS FOR SUCCESS

- Allow trial lenses to equilibrate for at least 10 minutes before assessing fit and vision
- Confirm axis orientation
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks



REFINE IF NEEDED: CONFIRM AXIS ORIENTATION ▶ DETERMINE EYE DOMINANCE ▶ FOLLOW GUIDANCE BELOW

	NEAR VISION			DISTANCE VISION			
		DOMINANT EYE	NON-DOMINANT EYE		DOMINANT EYE	NON-DOMINANT EYE	
IF PATIENT IS WEARING:	TWO LOW ADDS	Initial Lens	Low Add	Low Add	Initial Lens	Low Add	Low Add
		Refinement 1	Low Add	High Add	Refinement 1	Bausch + Lomb ULTRA® for Astigmatism	Low Add
		Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.			Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Bausch + Lomb ULTRA® for Astigmatism) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		
IF PATIENT IS WEARING:	TWO HIGH ADDS	Initial Lens	High Add	High Add	Initial Lens	High Add	High Add
		Refinement 1	High Add	Add +0.25D to the non-dominant eye	Refinement 1	Low Add	High Add
		Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.			Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		

BRING INNOVATION TO MORE PATIENTS

LENS PARAMETERS

MATERIAL:	samfilcon A	SPHERICAL POWERS:	+4.00D to -6.00D in 0.25D steps	
LENS MATERIAL TECHNOLOGY:	MoistureSeal® technology	CYLINDER POWERS:	-0.75D, -1.25D, -1.75D	-2.25D, -2.75D
WATER CONTENT:	46%	AXES:	10° to 180° in 10° steps	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°
OXYGEN PERMEABILITY:	114 Dk	ADD POWERS:	Low: +0.75D to +1.50D High: +1.75D to +2.50D	
LENS DESIGN TECHNOLOGY:	3-Zone Progressive™ Design, OpticAlign® Design	VISIBILITY TINT:	Light blue	
BASE CURVE:	8.6 mm	REPLACEMENT SCHEDULE:	Monthly	
DIAMETER:	14.5 mm			