

FITTING GUIDE

Biotrue® ONEday for PRESBYOPIA



TAKE A CLOSER LOOK

9 out of 10 ECPs agreed that a fitting guide made this 3-Zone Progressive™ Design easy to fit.¹ In fact, **80%** of patients were successfully fit in **one visit** and **96%** were successfully fit in **two visits**.*

*When the ECP followed the fitting guide for PureVision®2 for Presbyopia contact lenses. Thirty-nine ECPs (from 10 countries) refitted 441 existing soft contact lens wearing presbyopes into PureVision®2 Presbyopia lenses. Patients returned for follow-up visits after 1-2 weeks. ECP assessment of lens performance including ease of fit, and patient satisfaction with lenses in real-world conditions, were measured using a 6-point agreement survey.

BAUSCH+LOMB
See better. Live better.

3-ZONE PROGRESSIVE™ DESIGN

Seamless transitions between the three zones

NEAR


91% of patients report clear vision²

INTERMEDIATE

96% of patients report clear vision²

DISTANCE

95% of patients report clear vision²



2 ADD POWERS

Simplified fit with a low Add up to +1.50D and a high Add up to +2.50D

PATENTED DEHYDRATION BARRIER

Maintains nearly 100% of lens moisture for a full 16 hours¹

ADD SELECTION: FOR INITIAL LENS FITTING

SPECTACLE ADD	BOTH EYES
+0.75D to +1.50D	Low Add
+1.75D to +2.50D	High Add

SUGGESTED PATIENT CRITERIA

- Good motivation and realistic expectations
- Refractive astigmatism no greater than -1.00D

Select initial lenses

- Update spectacle refraction and Add power
- Determine ocular dominance for distance vision
- Select contact lens distance prescription based on spherical equivalent from spectacle Rx, adjusted for vertex distance if necessary
- Choose trial contact lenses based on the above calculation and select Add

Evaluate initial lenses

- Allow trial lenses to equilibrate for at least 10 minutes before assessing fit and vision
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks

NEAR VISION

IF PATIENT IS WEARING:

	DOMINANT EYE		NON-DOMINANT EYE	
TWO LOW ADDS	Initial Lens	Low Add	Low Add	
	Refinement 1	Low Add	High Add	
	<p>Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			
	DOMINANT EYE		NON-DOMINANT EYE	
TWO HIGH ADDS	Initial Lens	High Add	High Add	
	Refinement 1	High Add	Add +0.25D to the non-dominant eye	
	<p>Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using handheld lenses, and continue evaluating vision binocularly at normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			

DISTANCE VISION

IF PATIENT IS WEARING:

	DOMINANT EYE		NON-DOMINANT EYE	
TWO LOW ADDS	Initial Lens	Low Add	Low Add	
	Refinement 1	Biotrue® ONEday sphere	Low Add	
	<p>Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Biotrue® ONEday spherical lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			
	DOMINANT EYE		NON-DOMINANT EYE	
TWO HIGH ADDS	Initial Lens	High Add	High Add	
	Refinement 1	Low Add	High Add	
	<p>Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			

Fast, predictable, fitting with just 2 ADDs

	+0.75D	+1.00D	+1.25D	+1.50D	+1.75D	+2.00D	+2.25D	+2.50D
Biotrue® ONEday for Presbyopia	LOW ADD				HIGH ADD			
Air Optix Aqua Multifocal	LOW ADD			MED ADD		HIGH ADD		
1-Day Acuvue Moist Multifocal	LOW ADD			MID ADD		HIGH ADD		
Biofinity Multifocal (D/N)	+1.00D		+1.50D		+2.00D		+2.50D	

Lens parameters

Biotrue® ONEday for Presbyopia	
MATERIAL:	nesofilcon A
LENS MATERIAL TECHNOLOGY:	Surface Active Technology™
WATER CONTENT:	78%
OXYGEN TRANSMISSION:	42 Dk/t @ -3.00D
LENS DESIGN TECHNOLOGY:	3-Zone Progressive™ Design, center-near aspheric optics
BASE CURVE:	8.6 mm
DIAMETER:	14.2 mm
CENTER THICKNESS:	0.1 mm @ -3.00D
SPHERICAL POWERS:	+6.00D to -9.00D (in 0.25D steps, including plano)
ADD POWERS:	Low: +0.75D to +1.50D spectacle Add High: +1.75D to +2.50D spectacle Add
VISIBILITY TINT:	Light blue
MODALITY:	Daily disposable; Daily wear indication
UVA/UVB*:	Yes

***WARNING:** UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. The effectiveness of wearing UV-absorbing contact lenses in preventing or reducing the incidence of ocular disorders associated with exposure to UV light has not been established at this time. You should continue to use UV-absorbing eyewear as directed. NOTE: Long-term exposure to UV radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors, such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of outdoor activities). UV-blocking contact lenses help provide protection against harmful UV radiation. However, clinical studies have not been done to demonstrate that wearing UV-blocking contact lenses reduces the risk of developing cataracts or other eye disorders.

REFERENCES: 1. Data on file. Bausch & Lomb Incorporated. Rochester, NY. 2. Results of an online survey with patients who participated in a trial experience program with Biotrue® ONEday for Presbyopia contact lenses and wore their trial lenses for 4 or more days (n=261). Survey questions were top 3-box scores (% Strongly Agree, Agree, Slightly Agree) on a 6-point agreement scale, with a margin of error +/- 3.4%.

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. All other products/brand names and/or logos are trademarks of the respective owners.
©2019 Bausch & Lomb Incorporated or its affiliates. MTB.0059.USA.19

BAUSCH + LOMB
See better. Live better.